# PROTHESE D'EPAULE SHOULDER PROSTHESIS

# NOTICE D'INSTRUCTIONS INSTRUCTIONS FOR USE

### FR. FH ORTHOPEDICS

3 rue de la Forêt 68990 HEIMSBRUNN - FRANCE

Tél. +33 3 89 81 90 92 / Fax: +33 3 89 81 80 11 orthopedie@fhorthopedics.fr / www.fhorthopedics.fr

Consult instructions for use

R

### USA. FH ORTHO INC.

4908 N. Elston, Chicago, Illinois 60630 - USA Tel.: +1 773 290 1039 / 844-77 FHINC

Fax: +1 (773) 539 9328

info-us@fhorthopedics.com / www.FHortho.com

### UK. FH ORTHO LTD

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# Instructions for use - SHOULDER PROSTHESIS

# WARNING AND INDICATIONS

#### IT IS IMPORTANT TO CAREFULLY READ THIS LEAFLET

If you need further information, do not hesitate to contact the distributor's sales department. This document is not exhaustive. Neither is it an operative technique manual giving implantation details. In all cases, refer to the documents supplied

en

Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

FH ORTHOPEDICS supplies a range of shoulder implants:

- simple humeral prosthesis: humeral stems (press-fit or cemented) and humeral heads (centred
- total anatomical prosthesis (cemented glenoid implant with 4 pegs); humeral stems (press-fit or cemented), humeral heads (centred and off-centred) and cemented glenoids.
- total anatomical prosthesis (porous glenoid implant): humeral stems (press-fit or cemented). humeral heads (centred and off-centred), glenoid insert and porous glenoid base (cementless)
- total reverse prosthesis: humeral stems (press-fit or cemented), STD humeral insert, glenosphere and metal-back or porous glenoid base (without cement)
- All the implants of shoulder prosthesis are used for primary or revision surgeries.

The humeral stem Ø6 is intended for cemented application. The other humeral stems are intended for cemented or cementless application.

The glenoid bases (metal-back or porous) are intended for cementless application with the addition of cortical or cancellous bone screws. The porous glenoid base must be used for a total anatomical or reverse prosthesis. The metal-back glenoid base must be used only for reverse prosthesis.

#### Compatibility for prosthesis:

For anatomical shoulder prosthesis, the ideal mismatch (a radial mismatch of 4mm) between the head and the glenoid is obtained with same size of implants (same colour on the labelling). Ex: Ø46 with alenoid 46

For a cementless glenoid, the glenoid insert must have the same size as the glenoid base

Glenoid insert	Porous glenoid base
44	44S / 44S-LP / 44 / 44-LP
46	46 / 46-LP
48	48

For reverse shoulder prosthesis, the glenosphere and the humeral insert must have the same diameter.

le losphere and gleriold base must be assembled as follows.					
glenosphere	Metal-back glenoid base	Porous glenoid base			
Ø36	44 / 44S / 44R	44S / 44S-LP / 44 / 44-LP			
Ø39	44 / 44S / 44R / 46	44S / 44S-LP / 44 / 44-LP / 46 / 46-LP			
Ø42	46 / 48	46 / 46-LP / 48			

Cortical and cancellous bone screws must have the same material as the implant which will be used.

#### IMPLANT MATERIALS:

The implant materials are mentioned on each product label

Stainless steel and cobalt-chromium alloy (CoCr) usually contains nickel.

The implant is radio-opaque except if made out of UHMWPE.

### PERFORMANCES:

The device enables to restore the functional articular anatomy and the balance between the soft tissues. It ensures stable fixation and bone-saving stock.

Cortical and cancellous bone screws ensure fixation and holding of the implant on the bone and consolidation of the fracture.

#### INDICATIONS FOR USE:

The prostheses from FH Industrie are designed for specific indications such as:

- SIMPLE HUMERAL PROSTHESIS
- Fracture dislocation or complex four part fracture of the proximal humerus - Humeral head necrosis without injury to the glenoid cavity.
- Extensive humeral head cartilage damage without injury to the glenoid cavity
- Centred osteoarthritis with a glenoid cavity not allowing implantation of a glenoid implant.
- Rheumatoid polyarthritis with thin rotator cuff.
- Off-centred osteoarthritis with irreparable cuff, and with maintained active elevation of at least

TOTAL ANATOMICAL PROSTHESIS (CEMENTED GLENOID IMPLANT WITH 4 PEGS)

- Centred glenohumeral osteoarthritis with functional rotator cuff
- Rheumatoid polyarthritis with functional rotator cuff
- Fracture seguela, functional rotator cuff with glenoid injury

TOTAL ANATOMICAL PROSTHESIS (POROUS GLENOID IMPLANT)

- Centred glenohumeral osteoarthritis
- Rheumatoid polyarthritis
- Post-traumatic sequela with glenoid injury
- Fractures of the proximal humerus with glenoid injury
- Revision for glenoid loosening
- Glenoid bone loss, where bone graft is needed
- A functional rotator cuff is necessary to use this device

The porous glenoid base is intended for cementless application with the addition of bone screws for fixation.

REVERSE PROSTHESIS (METAL-BACK OR POROUS GLENOID IMPLANT)

The ARROW Reverse Shoulder Prosthesis is indicated for patients with severe shoulder arthropathy and a grossly deficient rotator cuff or a previously failed shoulder joint replacement with a grossly deficient rotator cuff. A functional deltoid muscle and adequate glenoid bone stock are necessary to use this device. The humeral stem is intended for cemented or cementless application while the glenoid baseplate (metal-back or porous) is intended for cementless application with the addition of bone screws for fixation.

# CONTRAINDICATIONS:

- Infection or latent infection
- Mental or neuromuscular disorders that may constitute an unacceptable risk to the patient and be a source of postoperative complications.
- Insufficient bone stock
- Obesity that could subject the prosthesis to stresses which may adversely affect device fixation or the device itself
- Known allergy to one of the components of the material mentioned on the product label
- Metabolic diseases that may compromise bone regrowth
- Drug-dependency

The humeral stem Ø6 must not be used if the patient body weight exceeds 76kg.

### PRECAUTIONS:

A - GENERAL PRECAUTIONS

- The lifetime of the implants depend on many biological, biomechanical and other factors. Therefore, strict compliance with the indications, contraindications, precautions and warnings is essential.
- The results of arthroplasty depend on the orthopaedic history of the patient. The psychological preparation of the patient is essential.
- The patients must be informed of the limitations of the device, including the impact of excessive stresses due to weight and intense activity. The patients must be given advice on how to adapt their activity. The prosthesis cannot restore the functions of a normal healthy bone structure. The patient should not nurture unrealistic hopes as to functionality and must consult the surgeon in case of abnormal disorders in the operative area.
- Before clinical use, the surgeon and operating theatre personnel must familiarize with the device and associated ancillary material. They must have an understanding of all the aspects of the surgical operation and the limitations of the device.
- It is important to protect the components and all polished or coated bearing surfaces from abrasion, scratches or any other adverse effect of metal and abrasive objects.
- The bearing surfaces must be clean and free from pieces of bone, cement, etc., before implantation. The surgeon should make sure that the implant is in full contact with the bone, and should remove any excess cement

- Make sure that the presence of other devices will not impair the functional integrity of the
- Do not mix implants of various sources. The only possible component combinations are indicated above (§ Compatibility for prosthesis).

### B - CEMENTLESS IMPLANTS

During handling, the utmost care should be paid to the coating. Avoid contact with the coating before implantation, Biologically fixed implants (identification of product level; HAP and/or T40) must not be cemented

## IMPLANTATION TECHNIQUE:

Specific instrumentation, delivered nonsterile, is necessary for the placement of implants designed by FH Industrie and distributed by FH ORTHOPEDICS or their distributors. The operative techniques are available from FH ORTHOPEDICS or their distributors. For some implants. FH ORTHOPEDICS or their distributors can supply X-ray templates for the surgeon to select the implant size. The centred drilling of the glenoid base 44 R must be controlled with X-ray imaging.

Patients receiving an implant must be informed that the lifetime of the implant may depend on their weight and activity.

- Damaged or defective implants must not be used,
- The reuse of the device is forbidden because of its conception and can modify its mechanical characteristics.
- Using an implant of inappropriate size may reduce resistance to strain.
- The components belong to a modular system. FH Industrie will accept no responsibility if components from another manufacturer are used or if their components are put to other than their intended use
- The implants must not be modified or subjected to any treatment.
- If the centred bone loss of the glenoid does not allow an anchorage of at least 2/3 of the glenoid keel length in bone, a glenoid base 44 R or -LP must be used
- The size 44 S glenoid base can be used only with a good quality bone stock. The antero-posterior screw must be used with a poor press-fit.

<u>UNDESIRABLE SIDE EFFECTS:</u>
It is the surgeon's responsibility to provide the patient with all the information before the operation, and especially to inform the patient of:

- the risk of rupture of the implant following inappropriate activity, trauma, or other stresses specific to the patient's activity.
- the risk of the implant coming loose following insufficient initial fixation, latent infection, premature or excessive stressing, component mispositionning, or trauma,
- the risk of allergy to one of the components of the material mentioned on the product label.
- side effects that may require further operation, or revision,
- possible progressive and sometimes asymptomatic bone resorption which may occur around implant due to reaction to foreign bodies.

IMPORTANT INFORMATION FOR THE SURGEON: Bone resorption may occur as a natural consequence of device implantation. Wear particles, especially from UHMWPE components, may cause osteolysis that may require re-operation.

The most common side effects of an arthroplasty are the following:

- Prosthesis dislocation
- Haematoma, late healing, deep vein thrombosis, pulmonary thrombosis, and blood vessel injury
- Pulmonary embolism
- Infection

with a metallic implant

- Cardiovascular problems

This device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of this device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. Before an MR exam, the patient should report to professional healthcare that he has been implanted

# PACKAGING AND STERILIZATION

The product label information gives traceability of manufacture. All implants are delivered sterile. Check for perfect sealing of packaging items (peel pouches or shells and seals) and overall integrity before using the implants. Do not use a product with a damaged package or a broken tamper-proof label. The resterilization of the device is forbidden (with any method) and can modify its mechanical

The stick-on dot on the outer package, that ensures the sterilization, must be red for radiation sterilization, or green for sterilization by ethylene oxide. This colour may be altered by adverse storage conditions: heat, humidity, light, etc.

In all cases, a stick-on dot that is yellow or purple (colour before sterilization) may indicate an unsterile product. In this case, the product must be returned. Desterilized implants are not taken back. The useby date is indicated on the product label.

# **DEVICE DISPOSAL**

When explanted, the implant must be handed over to a specialized service for environmentally safe disposal in compliance with strict hygiene rules. A product explanted because of a defect must be returned to the manufacturer after decontamination.