VIS BEPOD	NOTICE D'INSTRUCTIONS
BEPOD SCREW	INSTRUCTIONS FOR USE
BEPOD SCHRAUBEN	GEBRAUCHSANWEISUNG
TORNILLO BEPOD	INSTRUCCIONES DE USO
VITE BEPOD	ISTRUZIONI PER L'USO
PARAFUSO BEPOD	FOLHETO DE INSTRUÇÕES
ŚRUBY BEPOD	INSTRUKCJA STOSOWANIA
ŞURUBURI BEPOD	MANUAL DE UTILIZARE
BEPOD براغي	تعليمات الاستعمال
BEPOD SCHROEF	GEBRUIKSAANWIJZING

GROUP
FH
ORTHO™

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en

Instructions for use BEPOD SCREW

WARNINGS AND INDICATIONS FOR USE PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical techniques which show details about implantation. In all cases, please consult the included documents.

DESTINATION

The screws in the bePOD range are intended for implantation in the foot.

DESCRIPTION

FH Orthopedics supplies the following bePOD screw ranges:

- 3A cannulated screw
- Percutaneous screw
- 45B bevelled cannulated screw
- Arthrodesis cannulated screw
- Scored self-compressive screw
- EZ Weil and EZ Weil Power drive screws
- Hindfoot cannulated screw Ø7.2
- Hindfoot cannulated screw Ø5

BePOD screw range

EZ WEIL SCREW / EZ WEIL POWER DRIVE SCREW			
Length (mm)	10	12	14
Diameter (mm)			
Ø2.5	✓	✓	✓

45B BEVELLED CANNULATED SCREWS																							
Length (mm)	1	1	2	2	2	2	2	3	3	3	3	3	4	4	4	4	4	5	5	5	5	6	6
Diameter (mm)	6	8	0	2	4	6	8	0	2	4	6	8	0	2	4	6	8	0	2	4	6	8	0
Ø3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

3A CANNULATED SCREW														
Length (mm)	10	12	14	16	18	20	22	24	26	28	30	32	33	34
Diameter (mm)														
Ø2.5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	×	×	×
Ø3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

PERCUTANEOUS CANNULATED SCREW												
Length (mm)	10	12	14	16	18	20	22	24	26	28	30	32
Diameter (mm)												
Ø2.5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	×
Ø3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

HINDFOOT CANNULATED SCREW Ø7.2														
Length (mm)	45	50	55	60	65	70	75	80	85	90	95	100	105	110
	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

HINDFOOT CANNULATED SCREW Ø5																	
Length (mm)	24	26	28	30	32	34	36	38	40	45	50	55	60	65	70	75	80
	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

ARTHRODESIS CANNULATED SCREW										
Length (mm)				28	30	32	34	36	38	40
Diameter (mm)										
Ø4				✓	✓	✓	✓	✓	✓	✓

SCORED SELF-COMPRESSIVE SCREW																								
Length (mm)											10	12	14	16	18	20	22	24	26	28	30	32	34	
Diameter (mm)																								
Ø3											✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are noted on each product label.

Stainless steel typically contains nickel.

The implant is radiopaque unless it is made of UHMWPE.

COMPATIBILITY BETWEEN DEVICES

The ancillary device to be used for insertion of the implants is specified in the surgical technique.

The hindfoot cannulated screws Ø5 can be used with the Calcanail® nail.

PERFORMANCE

These implants bring closer together and fix bone fragments in order to allow osteosynthesis.

INDICATIONS FOR USE:

3A cannulated screws and scored self-compressive screws are designed for the following indications:

- Treatment of Hallux Valgus of the first metatarsal shaft (big toe) by compressive osteosynthesis, following distal and proximal metatarsal and phalangeal osteotomy.
- Compressive osteosynthesis of fractures of several bones in the forefoot.

Percutaneous cannulated screws, 45B bevelled cannulated screws are designed for the following indications:

- Treatment of Hallux Valgus of the first metatarsal shaft (big toe) by non-compressive osteosynthesis, following distal and proximal metatarsal and phalangeal osteotomy.

Non-compressive osteosynthesis of fractures of several bones in the forefoot.

EZWeil and EZWeil Power Drive screws are designed for the following indications:

- Treatment of metatarsalgia on the lateral metatarsal shafts (side toes), by Weil osteosynthesis, following osteotomies carried out:
- by the opened surgery (standard Weil) technique,
- by the DMMO percutaneous technique (Distal Metatarsal Mini Invasive Osteotomy).

Hindfoot cannulated screws Ø7.2 and Ø5 are designed for the following indications:

- Osteosynthesis on the bones of the hindfoot
- Fixation of mono and bi-cortical osteotomies
- Stabilisation of various fractures (e.g.: calcaneal, talus, etc.)
- Subtalar and talocrural arthrodesis.

Arthrodesis cannulated screws are designed for the following indication:

- Arthrodesis of the metatarsal-phalangeal joint of the big toe.

CONTRAINDICATIONS

- 1 - Infection or latent infection: they are diagnosed in patients with one or several of the following symptoms:
 - A- Fever and/or local inflammation
 - B - Accelerated destruction of the joint surface or osteolysis visible on the x-ray
 - C- Unexplained high erythrocyte sedimentation rate
 - D- Elevated white blood cell count or changes observed during different patient monitoring tests.
- 2 - A mental or neuromuscular disorder that might create an unacceptable risk of instability, prosthetic fixation failure, or complications after surgery.
- 3 - Bone weakened by illness or infection and unable to provide sufficient support and fixation. Skeletal immaturity.
- 4 - Obesity that might put stress on the prosthesis and comprise either fixation of the device or the device itself.
- 5 - Known allergy to any any compound of the material mentioned on the product label.
- 6 - Overweightness.
- 7 - Metabolic diseases that might compromise bone regrowth.
- 8 - Drug addiction.

PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.
- The outcome of the osteosynthesis depends on the patient's medical history. It is indispensable that the patient be psychologically prepared.
- Patients must be informed of the limitations of osteosynthesis, including among others, the impact of overburdening caused by excessive weight or activities. They must be advised about how to modify their activities accordingly. Osteosynthesis cannot, under any circumstances, restore the functions previously provided by a normal, healthy segment. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.
- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.
- It is advisable to protect the components and all polished or coated surfaces against the risk of abrasion, scratching and all other harmful effects associated with metal or abrasive objects.
- The contact zones of the implant must be clean and free of bone debris before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.
- Check that the presence of other devices will not affect the integrity or operation of the device.
- Do not combine implants from various origins. The only possible combinations of components are indicated in the section COMPATIBILITY BETWEEN DEVICES.

IMPLANTATION TECHNIQUE

The specific ancillary equipment, delivered non sterile, is required for the implantation of screws designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FR ORTHOPEDICS or its distributors. It is important that only this ancillary equipment be used for implantation.

CONDITIONS OF USE

The implants must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS

- 1- During handling operations, the implant must be protected from impacts that may cause damage. Never use a part that is damaged, defective or has suffered an impact.
- 2- Exclusive use of the ancillary equipment provided by F.H.ORTHOPEDICS or its distributors is essential to ensure proper implantation of the implant.
- 3- The patient must be warned that intense activity and/or overweightness increases the strain on the implant.
- 4- It is essential that the screw is tightened with care: it must not be overtight (damage to the bone), or not tight enough (inadequate fixation).
- 5- The implant cannot be re-sterilised by any process.
- 6- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.
- 7- Infection sites or a potential infection away from the operation area may spread to the operation site and must be treated as a priority.
- 8 - The implants must not be modified or subjected to any treatment.
- 9- Non-compliance with the above recommendations relieves the manufacturer and distributor of all responsibility.
- 10- Pre-surgery X-rays must be examined to avoid any risk of collision with a previously implanted device.

ADVERSE EFFECTS

The surgeon is responsible for providing the patient with complete information before the surgical procedure, and in particular for informing them of the following risks:

- Implant breakage following inappropriate activities, trauma or other strain generated by the patient.
- Implant disassembly due to initially inadequate fixation, latent infection, premature or excessive loading, incorrect component positioning or trauma.
- Risk of allergy to one of the components of the material mentioned on the product label (adverse effects may require a second operation or revision).
- Potential for progressive asymptomatic bone resorption around an implant due to a foreign body reaction

These adverse effects may lead to another operation or a revision.

IMPORTANT INFORMATION TO GIVE THE SURGEON AND THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

During MRI exams, the patient must inform the operator that he/she has a metallic implant:

- Metallic implants carry a risk of migration and/or heating during an MRI exam.
- Metallic implants may cause artifacts (distortion and/or loss of signal) which can affect interpretation of the imaging examinations.

PACKAGING AND STERILIZATION

The information on the product label ensures the traceability of its manufacturing.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly closed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be resterilized by any method whatsoever at the risk of altering its mechanical characteristics. The sterilization indicator sticker on the outer packaging should be red for radiation sterilization (R) or green for EtO gas sterilization (EO). This colour may be altered by poor storage conditions: heat, humidity, light etc. In all cases, an indicator showing one of the pre-sterilization colours, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilized implants cannot be returned. The expiry date is noted in the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialized waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING

EZ Weil screw: 2011

EZ Weil Power Drive screw: 2016

45B bevel cannulated screw: 2016

3A cannulated screw: 2007

Percutaneous cannulated screw: 2007

Cannulated screw for the hindfoot Ø7.2: 2007

Cannulated screw for the hindfoot Ø5 : 2011

Arthrodesis cannulated screw: 2005

Scored self-compressive screw: 2001

The reference text is the French text.

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يمكن استعمال البراعي المعجوة لمؤخرة الرصغ بقطر 5 بالترافق مع البرغي ©Calcanail.
UK, FH ORTHO LTD