GROUP  
**FH**  
ORTHO™

**FH** ORTHOPEDICS

**PROTHESE DISCALE LP-ESP**

**LP-ESP DISC PROSTHESIS**

**BANDSCHEIBENPROTHESE LP-ESP**

**PRÓTESIS DISCAL LP-ESP**

**PROTESI DISCALE LP-ESP**

**PRÓTESE DISCAL LP-ESP**

**LP-ESP DISCUSPROTHESE**

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**WARNINGS AND INDICATIONS FOR USE**

**PLEASE READ THIS DOCUMENT CAREFULLY**

**Do not hesitate to consult your supplier's sales department if you need further information. This document is not exhaustive, nor is it a technical user manual containing detailed layout drawings. In any case, refer to the documents provided.**

**DESTINATION:**

The devices covered in these instructions for use are intended to be implanted at the level of the lumbar spine between L2 and S1 during lumbar disc arthroplasty.

**DESCRIPTION:**

FH ORTHOPEDICS supplies the LP-ESP lumbar disc prosthesis. It is composed of two metal end plates which are separated by a deformable portion that is intended to allow the disc to move.

**MATERIALS THE IMPLANTS ARE MADE FROM**

The materials that the implants are made from are noted on each product label.

The implant is radio-visible, except for the polycarbonate-urethane and the silicone.

**COMPATIBILITY BETWEEN DEVICES:**

This is a monoblock device and does not require additional fixation elements.

**PERFORMANCE:**

The LP-ESP disc prosthesis must enable pain reduction, re-establish the lumbar curvature, and reproduce the functions of the disc. This medical device is used for primary treatment.

**INDICATIONS FOR USE:**

The FH ORTHOPEDICS LP-ESP lumbar disc prosthesis is designed to address specific indications described in the marketing documents, such as:

- lower back pain caused by lumbar disc disease that is resistant to medical treatment
- lower back pain caused by disc disease after treating a herniated disc
- lumbar radiculopathy caused by a herniated disc recurrence (apart from a sequestered disc)

**CONTRAINDICATIONS:**

- 1- Spinal stenosis, radiculopathy
- 2- Significant segmental instability
- 3- Spinal deformation, spondylolisthesis greater than 25%
- 4- X-ray confirmation of severe lesions or degeneration of the facet joints

- 5- Osteoporosis, osteochondrosis, and severe osteopenia
- 6- Chronic or severe, local, spinal, or systemic infections
- 7- Metabolic and systemic diseases
- 8- Pathologies and surgical situations that preclude any benefit of spinal surgery
- 9- Sensitivity to the implant's materials
- 10- Dependency with regards to medications: drug addiction or alcoholism
- 11- Pregnancy
- 12- Obesity
- 13- Lack of cooperation from the patient

## **PRECAUTIONS:**

### **A - GENERAL PRECAUTIONS**

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.

The outcome of intervertebral disc prosthesis depends on the patient's medical history. It is indispensable that the patient be psychologically prepared.

Patients must be informed of the limitations of the device, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.

- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.

- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.

- It is advisable to make sure that the presence of other devices does not compromise the device's functional integrity.

- Implants provided by FH ORTHOPEDICS must not be mixed with implants of another origin, or with different systems.

### **B - CEMENTLESS IMPLANTS**

During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants (identification on the product label: T40 + HAP) must not be cemented.

### **IMPLANTATION TECHNIQUE:**

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FH ORTHOPEDICS or its distributors. For certain implants, FH ORTHOPEDICS or its distributors may provide x-ray templates to enable the surgeon to evaluate the size of the implant.

### **CONDITIONS OF USE:**

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. The implantation must be performed with suitable non-damaged instrumentation by applying the recommendations in the surgical technique and instructions for use, to treat a patient with the indications defined above, either as a primary or revision arthroplasty procedure.

### **WARNINGS:**

- Patients receiving an implant must be informed that the implant's lifespan is affected by their weight and activity level
- Damaged or defective implants must not be used.
- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.
- - Using an implant that is of the wrong size or sacral tilt runs the risk of decreasing its stress resistance.
- FH Industrie is not responsible for use of another manufacturer's components, or use of components for purposes other than those intended.
- The implants must not be modified or subjected to any treatment.

### **ADVERSE EFFECTS:**

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

- Implant breakage following inappropriate activities, trauma or other loads related to the patient's activities
- Implant disassembly due to initially inadequate fixation, latent infection, premature or excessive loading, incorrect component positioning or trauma
- Allergy to one of the compounds of the material listed on the product label
- Side effects that could require a new surgical procedure or revision
- Potential for progressive asymptomatic bone resorption around prosthetic components due to a foreign body reaction
- Risks of degradation of adjacent vertebrae and risk of secondary movements

### **IMPORTANT INFORMATION FOR THE SURGEON:**

**Any serious accident that occurs in relation to the device should be**

**reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.**

If undergoing an MRI exam, the patient must inform the staff that he/she has a metal implant:

- Metal implants carry a risk of migration and/or heating during an MRI exam
- Metal implants can cause artifacts (signal distortion and/or signal loss) which may interfere with the interpretation of the imaging results.

**PACKAGING AND STERILISATION**

The information on the product label ensures manufacturing traceability and indicates which sterilization method has been used (if applicable).

**IMPLANTS DELIVERED STERILE:**

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly sealed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. The product must not be resterilized by any method whatsoever at the risk of altering its mechanical characteristics.

The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation. In all cases, an indicator showing the pre-sterilisation colours, yellow, may indicate that the product is not sterile and, in this case, the product must be returned.

Desterilised implants cannot be returned.

The expiry date is noted in the product label.

**DEVICE DISPOSAL**

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

**DATE OF FIRST CE MARKING:**

LP-ESP Disc prosthesis: 2007

**The reference text is the French text.**

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**UK, FH ORTHO LTD**