GROUP
FH
ORTHO

FH ORTHOPEDICS

**CUPULE BLINDEE PHARO
TETE CEPHALIQUE CENTREE
TETE CEPHALIQUE UNIPOLAIRE CONE 12/14**

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Instructions for use
PHARO BIPOLAR SHELL – CENTERED CEPHALIC HEAD – UNIPOLAR CEPHALIC HEAD TAPER 12/14

WARNINGS AND INDICATIONS FOR USE
PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical techniques which show details about implantation. In all cases, please consult the included documents.

DESCRIPTION

FH ORTHOPEDICS offers a range of PHARO BIPOLAR SHELLS, a range of CENTERED CEPHALIC HEADS and a range of UNIPOLAR CEPHALIC HEADS TAPER 12/14.

PHARO BIPOLAR SHELL

SIZE FOR PHARO BIPOLAR SHELL												
External diameter (mm)	41	43	45	47	49	51	53	55	57	59	61	63
Internal diameter (mm)												
Ø22.2	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ø28	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

CENTERED CEPHALIC HEAD and UNIPOLAR CEPHALIC HEAD TAPER 12/14

SIZE FOR CENTERED CEPHALIC HEAD and UNIPOLAR CEPHALIC HEAD TAPER 12/14																									
Diameter (mm)	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61
CENTERED CEPHALIC HEAD	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓	x	✓
UNIPOLAR CEPHALIC HEAD TAPER 12/14	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Necklength (mm)	-5	-3.5	+0	+3.5	+7	+10
INTERMEDIATE CENTERED RING TAPER 12/14	✓	✓	✓	✓	✓	✓

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are noted on each product label.

Stainless steel typically contains nickel.

The implant is radiopaque unless it is made of UHMWPE.

COMPATIBILITY BETWEEN DEVICES

It is imperative to assemble:

- the « CENTERED CEPHALIC HEAD » with the « INTERMEDIATE CENTERED RING TAPER 12/14 ».
- the « INTERMEDIATE CENTERED RING TAPER 12/14 », and the « UNIPOLAR CEPHALIC HEAD TAPER 12/14 » with the femoral stems from FH INDUSTRIE SC (cementless) or AC (cemented) with a 12/14-definition TAPER.
- the « PHARO BIPOLAR SHELL » with heads of the same diameter (for example: Ø28).

Material compatibility of PHARO BIPOLAR SHELLS:

	PHARO BIPOLAR SHELL
Stainless steel heads TAPER 12/14 FH INDUSTRIE	✓
CoCr heads TAPER 12/14 FH INDUSTRIE	✓
Stainless steel skirted heads TAPER 12/14 FH INDUSTRIE	✓
CoCr skirted heads TAPER 12/14 FH INDUSTRIE	✓
Alumina heads BIOLOX® FORTE TAPER 12/14 FH INDUSTRIE	✓
Alumina heads BIOLOX® DELTA TAPER 12/14 CERAMTEC	✓ see table below

FEMORAL HEADS MANUFACTURED BY CERAMTEC

Necklength (mm)			
BIOLOX® DELTA Ø28 TAPER 12/14	28 12/14 S Ref. 38.49.7175.445.00	28 12/14 M Ref. 38.49.7175.455.00	28 12/14 L Ref. 38.49.7175.465.00

The ancillary device to be used for the insertion of hip implants is specified in the surgical technique.

PERFORMANCE

The device helps restore joint function and stands the usual stresses to which the joint in question is submitted.

INDICATIONS FOR USE:

The range PHARO BIPOLAR SHELL, the range CENTERED CEPHALIC HEAD and the range UNIPOLAR CEPHALIC HEAD TAPER 12-14 are designed to meet the specific indications during intermediate hip replacement such as:

- Traumatic injuries
- Femoral neck fracture
- Osteonecrosis of the head

CONTRAINDICATIONS

1 - Infection, or latent infection: they are diagnosed in patients with one or several of the following symptoms:

- A- Fever and/or local inflammation
- B- Accelerated destruction of the joint surface or osteolysis visible on the x-ray
- C- Unexplained high erythrocyte sedimentation rate
- D- Elevated white blood cell count or changes observed during different patient monitoring tests.

2- A mental or neuromuscular disorder that might create an unacceptable risk of instability, prosthetic fixation failure, or complications after surgery.

3 - Insufficient bone stock.

4 - Known allergy to any component of the material mentioned on the product label.

5 - Metabolic diseases that might compromise bone regrowth.

6 - Drug addiction.

PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.

The outcome of arthroplasty depends on the patient's medical history. It is indispensable that the patient be psychologically prepared.

Patients must be informed of the limitations of the prosthesis, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised about how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.

- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.

- The contact zones of the implant must be clean and free of debris (bone) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.

- Care should be taken to ensure that the presence of other devices will not harm the functional integrity of the device.

- Do not combine implants from various origins. The only possible combinations of components are indicated in the section « Compatibility between devices ».

IMPLANTATION TECHNIQUE

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FR ORTHOPEDICS or its distributors. For certain implants, FH ORTHOPEDICS or its distributors may provide x-ray templates to enable the surgeon to evaluate the size of the implant.

CONDITIONS OF USE

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS

Patients receiving a prosthesis must be informed that the lifespan of the implant may depend on their weight and their activity level.

- Damaged or defective implants must not be used.
- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.
- Using an implant that is of the wrong size or lateralization runs the risk of decreasing its stress resistance.
- The components are part of a modular system. FH Industrie refuses any responsibility in the event that components from another manufacturer are used, or that its components are not used as they are intended to be.
- The implants must not be modified or subjected to any treatment.

ADVERSE EFFECTS

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

- pain
- loss of mobility
- wear of the UHMWPE components. Wear particles, especially from components made of UHMWPE, may cause osteolysis that may require another surgery.
- implant rupture as a result of inappropriate activities, trauma, or other stress specific to the patient's activities.
- loosening or dislocation of an implant following inadequate initial fixation, latent infection, early or excessive stress, poor positioning of the components, or trauma.
- allergy to one of the compounds of the material noted on the product label.
- possible evolving, sometimes asymptomatic, bone resorption that may appear around prosthetic components as a result of reactions to a foreign body.
- hematoma, slowed healing, deep vein thrombosis, pulmonary thrombosis, vascular lesions.
- pulmonary embolism
- infection, and possibly death
- a cardiovascular problem

These adverse effects may lead to another operation or a revision.

IMPORTANT INFORMATION TO GIVE THE SURGEON AND THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

Metallic implants carry a risk of migration and/or heating during an MRI exam. During MRI exams, the patient must inform the medical staff that he/she has a metallic implant.

PACKAGING AND STERILIZATION

The information on the product label ensures the traceability of its manufacturing.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly closed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be resterilized by any method whatsoever at the risk of altering its mechanical characteristics. The sterilization indicator sticker on the outer packaging should be red for radiation sterilization (R) or green for ETO gas sterilization (EO). This colour may be altered by poor storage conditions: heat, humidity, light, etc. In all cases, an indicator showing one of the pre-sterilization colours, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilized implants cannot be returned. The expiry date is noted on the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialized waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE FOR FIRST CE MARKING

PHARO BIPOLAR SHELL: 2006

Centered cephalic head: 2004

UNIPOLAR CEPHALIC HEAD TAPER 12/14: 2004

Intermediate centered ring TAPER 12/14: 2004

The reference text is the French text.

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