**NOTICE D'INSTRUCTIONS**  
**INSTRUCTIONS FOR USE**  
**INSTRUCCIONES DE USO**

CE 0459

Imprimé en France  
EMNI 2369 001 ind. B  
Revision 2017-07

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en Instructions for use – STS IMPLANTS

**WARNING AND INDICATIONS**

**IT IS IMPORTANT TO CAREFULLY READ THIS LEAFLET**

If you need further information, do not hesitate to contact the distributor's sales department. This document is not exhaustive. Neither is it an operative technique manual giving implantation details. In all cases, refer to the documents supplied

**CAUTION**

Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

**DESCRIPTION:**

FH ORTHOPEDICS supplies a range of shoulder implants:

- textile assembly,
- cementation assembly,
- fixation system

All the implants of shoulder prosthesis are used for primary or revision surgeries.

**Compatibility for prosthesis:**

The textile assembly for fracture must be used only with Ø8 humeral stems with logotype "stem STS".

The cementation assembly must be used only with Ø8 humeral stems and Ø8 fracture stems.

The fixation system must be used only with fracture stems.

**IMPLANT MATERIALS:**

The implant materials are mentioned on each product label. Stainless steel and cobalt-chromium alloy (CoCr) usually contains nickel.

The implant is radio-opaque except if made out of UHMWPE or PET.

**PERFORMANCES:**

The device enables to restore the functional articular anatomy and the balance between the soft tissues. It ensures stable fixation and bone-saving stock.

**INDICATIONS FOR USE:**

The prostheses from FH Industrie are designed for specific indications such as:

- Fracture of the upper extremity of the humerus

**CONTRAINDICATIONS:**

- Infection or latent
- Mental or neuromuscular disorders that may constitute an unacceptable risk to the patient and be a source of postoperative complications.
- Insufficient bone stock
- Obesity that could subject the prosthesis to stresses which may adversely affect device fixation or the device itself.
- Known allergy to one of the components of the material mentioned on the product label.
- Metabolic diseases that may compromise bone regrowth
- Drug-dependency

**PRECAUTIONS:**

**GENERAL PRECAUTIONS**

- The lifetime of the implants depend on many biological, biomechanical and other factors. Therefore, strict compliance with the indications, contraindications, precautions and warnings is essential.
- The results of arthroplasty depend on the orthopaedic history of the patient. The psychological preparation of the patient is essential.
- The patients must be informed of the limitations of the device, including the impact of excessive stresses due to weight and intense activity. The patients must be given advice on how to adapt their activity. The prosthesis cannot restore the functions of a normal healthy bone structure. The patient should not nurture unrealistic hopes as to functionality and must consult the surgeon in case of abnormal disorders in the operative area.
- Before clinical use, the surgeon and operating theatre personnel must familiarize with the device and associated ancillary material. They must have an understanding of all the aspects of the surgical operation and the limitations of the device.
- It is important to protect the components and all polished or coated bearing surfaces from abrasion, scratches or any other adverse effect of metal and abrasive objects.
- The bearing surfaces must be clean and free from pieces of bone, cement, etc., before implantation. The surgeon should make sure that the implant is in full contact with the bone, and should remove any excess cement.
- Make sure that the presence of other devices will not impair the functional integrity of the device.
- Do not mix implants of various sources. The only possible component combinations are indicated above (§ Compatibility for prosthesis).

#### **IMPLANTATION TECHNIQUE:**

Specific instrumentation, delivered nonsterile, is necessary for the placement of implants designed by FH Industrie and distributed by FH ORTHOPEDICS or their distributors. The operative techniques are available from FH ORTHOPEDICS or their distributors. For some implants, FH ORTHOPEDICS or their distributors can supply X-ray templates for the surgeon to select the implant size.

#### **WARNINGS:**

- Patients receiving an implant must be informed that the lifetime of the implant may depend on their weight and activity.
- Damaged or defective implants must not be used,
  - The reuse of the device is forbidden because of its conception and can modify its mechanical characteristics.
  - Using an implant of inappropriate size may reduce resistance to strain.
  - The components belong to a modular system. FH Industrie will accept no responsibility if components from another manufacturer are used or if their components are put to other than their intended use.
  - The implants must not be modified or subjected to any treatment.

#### **UNDESIRABLE SIDE EFFECTS:**

**It is the surgeon's responsibility to provide the patient with all the information before the operation, and especially to inform the patient of:**

- the risk of rupture of the implant following inappropriate activity, trauma, or other stresses specific to the patient's activity,
- the risk of the implant coming loose following insufficient initial fixation, latent infection, premature or excessive stressing, component mispositioning, or trauma,
- the risk of allergy to one of the components of the material mentioned on the product label,
- side effects that may require further operation, or revision,
- possible progressive and sometimes asymptomatic bone resorption which may occur around implant due to reaction to foreign bodies.

**IMPORTANT INFORMATION FOR THE SURGEON:** Bone resorption may occur as a natural consequence of device implantation. Wear particles, especially from UHMWPE components, may cause osteolysis that may require re-operation.

The most common side effects of an arthroplasty are the following:

- Prosthesis dislocation
- Haematoma, late healing, deep vein thrombosis, pulmonary thrombosis, and blood vessel injury
- Pulmonary embolism
- Infection
- Cardiovascular problems

During MRI/CT scan examinations, ask the patient if he or she wears a metal implant. This device has not been evaluated for safety and compatibility. This device has not been tested for heating or migration. The specialist in charge of the MRI/CT scan examinations must hold on security of device and patient, the artefact risk and the information loss.

#### **PACKAGING AND STERILIZATION**

The product label information gives traceability of manufacture and the method of sterilisation (if used).

#### **IMPLANTS DELIVERED NONSTERILE:**

- Before implantation, the implants must be washed, decontaminated, and packed for sterilization by steam autoclaving according to a process validated by the user.
- Recommended sterilization: cleaning and decontamination, then physical inactivation by steam in autoclave at a temperature  $\geq 134^{\circ}\text{C}$  for at least 18 minutes.

#### **IMPLANTS DELIVERED STERILE:**

- Check for perfect sealing of packaging items (peel pouches or shells and seals) and overall integrity before using the implants. Do not use a product with a damaged package or a broken tamper-proof label. The resterilization of the device is forbidden (with any method) and can modify its mechanical characteristics. The stick-on dot on the outer package, that ensures the sterilization, must be red for radiation sterilization, or green for sterilization by ethylene oxide. This colour may be altered by adverse storage conditions: heat, humidity, light, etc. In all cases, a stick-on dot that is yellow or purple (colour before sterilization) may indicate an unsterile product. In this case, the product must be returned. Desterilized implants are not taken back. The use-by date is indicated on the product label.

#### **DEVICE DISPOSAL**

When explanted, the implant must be handed over to a specialized service for environmentally safe disposal in compliance with strict hygiene rules. A product explanted because of a defect must be returned to the manufacturer after decontamination.

The reference text is the French text.

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