GROUP

FH
ORTHO[®]

FH ORTHOPEDICS

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Instructions for use STERILE PACKAGED AND DELIVERED INSTRUMENT

WARNINGS AND INDICATIONS FOR USE PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical techniques which show details about implantation. In any case, refer to the documents provided.

Under federal law (United States), the device can be sold, distributed and used by a doctor or a doctor's prescription.

DESCRIPTION

FH ORTHOPEDICS offers a range of sterile packaged and delivered instruments. These include temporary invasive devices such as reamers and guide pins.

The following instruments are available:

- 267483: Sterile reamer Ø4.5 mm
- 267482: Sterile pins Kit
- 266636: CoLS Disposable instruments, sterile
- 266637: CoLS 2.4 mm disposable guide pin, sterile (Lg 250 mm)
- 266638: CoLS 4.5 mm disposable drill, sterile
- 266639: CoLS 1.8 mm disposable screw pin, sterile (Lg 220mm)
- 263883: Set of two pins, sterile

References 266636, 266637, 266638, 266639, and 263883 are specific to and only available for the USA.

References 267483 and 267482 cannot be sold in the USA.

PERFORMANCE

The instrumentation is intended to be used to correctly insert the implant or implant system that it is associated with, in accordance with the surgical technique.

MATERIALS

The materials that the sterile delivered instruments are made from are noted on each product label.

They are manufactured from materials that comply with ISO 16061.

The instruments (or a fragment) are metallic, and therefore may be located using imaging equipment.

COMPATIBILITY

The instruments provided with these instructions are intended to be used alone or in combination with the implant or other instruments; see the surgical techniques of the devices in question to find out about the various combinations that may be implemented. They are available from FH ORTHOPEDICS or its distributors.

INDICATIONS

The ancillary equipment has been developed to meet specific indications. See the indications for the associated implant or implant system in the corresponding instructions for use.

CONTRAINdications

See the general contraindications of the associated implant or implant system.

PRECAUTIONS

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated all aspects of the surgery, as well as the limitations of the device.
- It is advisable to make sure that the presence of other devices does not compromise the instrument's functional integrity.
- It is advisable to take care to protect the instruments, particularly the cutting edges and the connector cords. Their condition affects the safe functioning of the instrument.
- Operating suite staff is responsible for checking, prior to the procedure, that F.H. Industrie instruments are connected correctly to the operating suite equipment available.
- Some ancillary equipment that includes temporary invasive devices such as pins and reamers must be used with caution. Rotation may give rise to conflict with the soft tissues or areas in contact. This equipment, tested by FH ORTHOPEDICS, or its distributors, should be handled with care. Any damage to the equipment may give rise to potential malfunction.
- The use of some instruments to be fitted on motors (drill bits, reamers, pins, etc.) may give rise to potential overheating between the instruments and the bone. It is recommended to spray these instruments with physiological saline solution during use.

CONDITIONS OF USE

The instruments must be used in an operating room by an orthopaedic surgeon familiar with this type of procedure under aseptic conditions while using good hygiene practices. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined in the instructions for use linked to the implant.

WARNINGS

- Damaged or defective instruments must not be used.
- The design of these instruments does not allow for their reuse, at the risk of altering their mechanical characteristics.
- The instruments must not be modified or subjected to any treatment.
- No chemical products (ex. alcohol) should be applied to the instrument during clinical practice.
- The ancillary equipment and the instruments featured in these instructions have been designed to be used exclusively with implants commercialised by FH ORTHOPEDICS. They cannot be used in combination with implants from other manufacturers.

Not heeding this warning may lead to deterioration of the associated instruments or implants and may also compromise their lifespan. The surgeon and the operating room staff are responsible for complications that may arise from not heeding this warning and not following the conditions of use for the ancillary equipment presented in these instructions.

ADVERSE EFFECTS

The surgeon is responsible for providing the patient with complete information before the surgical procedure, and in particular for informing them of the following risks:

Incorrectly using the ancillary equipment may lead

- to tissue damage,
- to premature wear,
- to deterioration of the instruments

This may represent a danger for the user, the patient, or a third party

These adverse effects may lead to another operation.

IMPORTANT INFORMATION:

Any serious accident that occurs in relation to the instrument should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

PACKAGING AND STERILISATION

The information on the product label ensures manufacturing traceability.

INSTRUMENTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly sealed and that all elements are undamaged before using the instruments. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be resterilised by any method whatsoever at the risk of altering its mechanical characteristics. The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation (R). This colour may be altered by poor storage conditions: heat, humidity, light etc. In all cases, an indicator showing one of the pre-sterilisation colours, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilised instruments cannot be returned. The expiry date is noted in the product label.

DEVICE DISPOSAL

In the event of a disposal that carries a particular or unusual risk, the instrument must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws.

Date of first CE marking: 2015, references below:

- 267483: Sterile reamer Ø4.5 mm
- 267482: Sterile pins Kit

The reference text is the French text.

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