GROUP
FH
ORTHO

FH ORTHOPEDICS

Notice d'instructions
Instructions for use
Gebrauchsanweisung
Instrucciones de uso
Istruzioni per l'uso
Folheto de instruções
Gebruiksaanwijzing
Φυλλάδιο οδηγιών
Instrukcja stosowania
Naudojimo instrukcija
Manual de utilizare

INSTRUMENT LIVRE STERILE
INSTRUMENT DELIVERED STERILE
STERIL VERPACKTES INSTRUMENT
INSTRUMENTO SUMINISTRADO ESTÉRIL
STRUMENTO FORNITO STERILE
INSTRUMENTO FORNECIDO ESTÉRIL
STERIEL GELEVERD INSTRUMENT
ΕΡΓΑΛΕΙΟ ΠΟΥ ΠΑΡΑΔΙΔΕΤΑΙ ΣΤΕΙΡΟ
PRZYRZĄD DOSTARCZONY STERYLNY
INSTRUMENTAS TIEKIAMAS STERILUS
INSTRUMENT LIVRAT STERIL

Consult
instructions
for use

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en

Instructions for use – INSTRUMENT DELIVERED STERILE

WARNINGS AND INDICATIONS FOR USE

PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical technique instructions which show details about implantation. In any case, refer to the documents provided.

CAUTION

Under federal law (United States), the device can be sold, distributed and used by a doctor or a doctor's prescription

DESCRIPTION:

FH ORTHOPEDICS supplies a range of instruments that are packaged and delivered sterile.

These are temporary invasive devices such as reamers and guide wires.

The following instruments are available:

- 267114 and 268017: stopped drill bit Ø5
- 267115: drill bit Ø3.2
- 268016: extended cutting guide pins Ø3 L90
- 269132: cannulated tapered reamer
- 269138: threaded pin Ø3 L170
- 269148 and 269149: cannulated drill bit Ø5
- 269240: quick-release peripheral drill Ø5
- 269478: metal-back instruments set
- 269479: full poly instruments set
- 269481: quick-release tapered reamer

Only the threaded pin ref 269138 is CE.

COMPATIBILITY OF DEVICE:

The instruments supplied with this IFU are intended to be used alone or in combination with the implant or with other instruments. Refer to the surgical technique associated with the devices in question to determine which combinations can be used. This information is available from FH ORTHOPEDICS or its distributors.

MATERIALS:

The materials that the instruments are made from are listed on each product label. They are manufactured from materials that conform to ISO 16061. Since the instruments (or a portion of them) are made of metal, they can be located during imaging.

PERFORMANCES:

The instrumentation is intended to be used to correctly insert the implant or the associated implant system, in accordance with the surgical technique.

INDICATIONS FOR USE:

The instrumentation sets were developed for specific indications. Refer to the indications in the IFU included with the implant or the associated implant system.

CONTRAINdications:

Refer to the general contraindications for the implant or the associated implant system.

PRECAUTIONS:

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its instrumentation. They must have fully understood all aspects of the surgery, as well as the limitations of the device.
- It is advisable to make sure that the presence of other devices does not compromise the instrument's functional integrity.
- The instrument, particularly its cutting edges and connectors, must be protected from damage. The instrument's integrity is essential to its safe use.
- Operating suite staff is responsible for checking, prior to the procedure, that F.H. Industrie instruments are connected correctly to the operating suite equipment available.

- Instrumentation sets consisting of temporary invasive devices such as K-wires and reamers must be used with care. Any rotation could lead to contact with soft tissues or the sterile drapes.
- This equipment, tested by FH ORTHOPEDICS, or its distributors, should be handled with care. Any damage to the equipment may give rise to potential malfunction.
- Use of instruments that are assembled with a surgical handpiece (e.g. drill bits, reamers, probes, milling cutters, screw taps, etc.) can cause overheating at the instrument–bone interface. It is recommended to spray these instruments with physiological saline solution during use.

CONDITIONS OF USE:

The instruments must be used in an operating suite, in an aseptic environment and according to good hygiene practices, by an orthopedic surgeon who regularly does this type of implantation surgery. Implantation must be carried out using the appropriate instrumentation set, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications listed in the implant's IFU.

WARNINGS:

- Damaged or defective instruments must not be used.
- This instrument's design does not allow it to be reused because its mechanical properties could be altered.
- The instruments must not be modified or subjected to any type of treatment.
- Do not apply chemical products (e.g. alcohol) to the instrument during clinical practice.
- The accessories and instruments featured in this IFU were designed solely to be used with implants sold by FH ORTHOPEDICS. They cannot be used with implants from other manufacturers.

Not complying with this warning can lead to the instruments or associated implants being damaged, which could compromise their lifespan. The surgeon and operating suite staff are responsible for any complications that may result from not complying with this warning and with the instructions listed in this IFU.

ADVERSE EFFECTS:

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

Incorrect use of the instrumentation can lead to:

- Tissue damage
- Premature wear
- Instrument damage.

This can be dangerous to the user, patient or a third-party.

These adverse effects may require that a new surgical procedure be performed.

IMPORTANT INFORMATION:

Any serious instrument-related accident should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient resides.

PACKAGING AND STERILIZATION

The information on the product label ensures the traceability of its manufacturing.

Verify that the packaging elements (peelable pouches or blister packs and films) are completely closed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be resterilized by any method whatsoever at the risk of altering its mechanical characteristics. The sterilization indicator sticker on the outer packaging should be red following radiation sterilization (R). This color may be altered by poor storage conditions: heat, humidity, light, etc. In all cases, an indicator showing one of the pre-sterilization colors, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilized instruments cannot be returned. The expiry date is listed on the product label.

DEVICE DISPOSAL

In the event of a disposal presenting a special or unusual risk, the instrument must be sent to a specialized waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws.

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