GROUPE

FH

ORTHO™

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Références.

- AAMI TIR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

EN - FH instrument operating, cleaning and sterilization instructions

Do not hesitate to consult your supplier's sales department if you need further information. This document is not exhaustive, nor is it a technical user manual containing detailed layout drawings. In any case, refer to the documents provided.

General

- This leaflet is applicable to all reusable medical devices manufactured by FH Industrie and distributed by FH ORTHO and/or its distributors and should be read attentively. It replaces FH instrument leaflets published prior to 2010.
 - Both new and old instruments should be processed with care in accordance with these instructions before use.
 - Single-use devices should not be resterilised with the exception of implants supplied non-sterile. These instructions are not applicable to FH medical devices supplied sterile which should not be reused.
 - Instrument cases do not provide a sterile barrier: sterilization packaging should be used to keep the device sterile.
 - The instructions provided in this leaflet have been validated by FH. It is the hospital's responsibility to ensure that reprocessing is carried out using suitable equipment and materials and that the staff in charge of reprocessing has received proper training to obtain the desired result. The equipment and processes require regular validation and monitoring.
 - Ancillary equipment is supplied non-sterile in transport containers. A transfer sheet (specifying the non-sterile status of the ancillary device) is supplied with each ancillary device for incoming quality control and prior to sterilization.
 - FH devices should exclusively be used for FH prosthesis implantation and explantation.
 - Stainless steel and cobalt-chrome (CoCr) alloy usually contain nickel. Nickel, cobalt and chrome can cause allergic reactions.
 - FHK-CK trial baseplates (references 268873 à 268878) are made of cobalt-chrome (CoCr). This alloy contains cobalt which is classified as CMR (Carcinogenic, Mutagenic or Reprotoxic). Nevertheless, the evaluation of the use of the cobalt-chrome alloy in FH medical devices has demonstrated their safety for the population in general.
-

Instructions for use

- A surgical procedure is available from the FH sales department or the distributor to provide surgeons with additional information on the recommended method for use. It is essential to read this document.
- Ancillary devices are designed to operate at ambient temperature. It is recommended not to use these devices directly after steam sterilization and to wait for the return to ambient temperature. Operating suite staff is responsible for checking, prior to the procedure, that F.H. Industrie instruments are connected correctly to the operating suite equipment available.
- Some ancillary devices consist of temporary invasive devices such as drills, reamers, milling cutters, screw taps and probes. These instruments should be used with caution. Rotation may give rise to conflict with the soft tissues or areas in contact.
- This equipment, tested by FH ORTHO, or its distributors, should be handled with care. Any damage to the equipment may give rise to potential malfunction.
- The equipment should only be used for its intended purpose in surgical procedures.
- The use of some instruments to be fitted on motors (drill bits, reamers, probes, milling cutters, screw taps, etc.) may give rise to potential overheating between the instruments and the bone. It is recommended to spray these instruments with physiological saline solution during use.

FH recommends processing all devices in accordance with the manual cleaning or combined manual/automatic cleaning instructions contained in this leaflet.

Warnings

- **The following precautions should be observed** by hospital staff working with contaminated or potentially contaminated medical devices. Pointed or sharp devices should be handled with the utmost caution.
- **Do not use metal brushes or scouring pads** during manual cleaning procedures. This equipment is liable to damage the instrument surface and finish. Use soft nylon bristle brushes and cleaning brushes of suitable sizes for the devices to be processed.
- To facilitate the subsequent steps, **do not allow contaminated devices to dry before reprocessing.**
- Physiological saline solution and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and **must not** be used.
- **Do not use mineral oil or silicone lubricants** as they coat micro-organisms, prevent direct contact of the surface with the steam and are difficult to remove.
- **Automatic cleaning using only a washer/disinfector is not sufficient for orthopaedic instruments.** A thorough manual cleaning before automatic cleaning procedure is required.
- pH-neutral cleaning and enzyme agents are recommended and preferred for cleaning reusable FH devices.
- When cleaning aluminium devices, only pH-neutral detergents should be used. Contact with strong alkaline detergents or with solutions containing sodium hydroxide, iodine or chlorine should be avoided as the aluminium may undergo chemical corrosion and the device may be damaged. The user should always refer to and apply the instructions provided by the cleaning product manufacturer.

Note: Drill bits, pins, rasps and sharp instruments may be treated with alkaline detergents. They should be inspected thoroughly after treating to ensure that the sharp edges are not damaged.

Instructions

- Remove excess bodily fluids and tissues on the instruments with a disposable lint-free pad. Place the devices in a tray of distilled water or cover with a damp cloth.
- Instruments **must** be cleaned within 30 minutes following use to minimise the risk of drying prior to cleaning.
- Used instruments **must** be transported to the procurement department in sealed or covered containers to prevent any further risk of contamination.
- If applicable, instruments with multiple components should be dismantled for effective cleaning. Take care not to mislay any small screws and components.
- All cleaning agents should be prepared in accordance with the dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning agents. It is important to use the recommended temperatures for optimal cleaning agent performances.

Note: Fresh cleaning solutions must be prepared.

The instruments **must** be removed from the metal or polymer trays when applying manual and/or automatic cleaning procedures. Lids, cases and instrument trays should be cleaned separately, except for non-sterile single-use implants. Plates and screws may remain in the tray or trolley for reprocessing.

- **The use of hard water should be avoided.** Softened tap water may be used for the initial rinse. The final rinses should be carried out with purified water to remove deposits on the instruments. One or more of the following processes may be used for water purification purposes: ultrafiltration (UF), reverse osmosis (RO), de-ionisation or equivalent.
- Steam (moist heat) is the validated sterilization method for FH instruments.

Manual disinfection/cleaning procedure

Step 1 Steep the instruments completely in an enzyme solution (ALKAZYM 0.5%) and leave to soak for 20 minutes. Use a soft nylon bristle brush to scrub the device gently until any visible stains have been removed. Pay

particular attention to rough areas, cavities, contact surfaces, connectors and other areas that are difficult to clean. Actuate hinged parts to access areas subject to difficult access.

Step 2 Remove the device from the enzyme solution and rinse with tap water for at least 3 minutes. Rinse any cavities, holes and any other areas subject to difficult access thoroughly and effectively with a pressurised spray.

Step 3 Place the cleaning agents in the ultrasonic cleaner. Steep the device completely in the cleaning solution and sonicate for at least 10 minutes at 45–50 kHz.

Step 4 Rinse the device in purified water for at least 3 minutes or until the blood or visible stains on the device or in the rinsing water have been removed. Rinse any cavities, holes and any other areas subject to difficult access thoroughly and effectively.

Step 5 Repeat the above ultrasonic cleaning and rinsing steps.

Step 6 Remove any excess moisture on the instrument with a clean, absorbent, lint-free cloth.

Manual/automatic disinfection/cleaning procedure

Step 1 Steep the instruments completely in an enzyme solution (ALKAZYM 0.5%) and leave to soak for 10 minutes. Use a soft nylon bristle brush to scrub the device gently until any visible stains have been removed. Pay particular attention to rough areas, cavities, contact surfaces, connectors and other areas that are difficult to clean. Actuate hinged parts to access areas subject to difficult access.

Step 2 Remove the devices from the enzyme solution and rinse with tap water for at least 3 minutes. Rinse any cavities, holes and any other areas subject to difficult access thoroughly.

Step 3 Place the instruments in a suitable washer/disinfector basket and process via standard instrument washer/disinfector cycle in accordance with the washer/disinfector manufacturer's instructions:

- Wash for 10 minutes at 93 °C minimum with a suitable washer/disinfector detergent solution.
- Rinse with deionised water
- Dry.

FH Industrie recommends the use of an ISO 15883 compliant washer/disinfector for instrument cleaning and disinfection.

Inspection before sterilization

- The cleaning procedures should be repeated until instruments are "thoroughly cleaned" and "visually clean".
 - Inspect each device carefully to ensure that any visible contamination has been removed. In the presence of contamination, repeat the cleaning/disinfection procedure.
 - Check the action of moving parts (e.g., hinges, connectors, sliding parts, etc.) throughout their range.
 - Check that there is no instrument deformation (particularly instruments involving connection to motors).
 - When the instruments have been assembled, check that the devices are assembled correctly with the components in contact.
 - Inspect devices for wear and tear: the device is not suitable for reuse if it presents corrosion, discoloration, staining or damage. Contact your sales representative to replace such defective device.
-

Sterile packaging

Commercially available medical grade steam sterilization wrapping FDA-cleared or pouches may be used to package individual instruments. The wrapping should be prepared according to the DIN EN ISO/ANSI AAMI ISO 11607-1/2 protective double wrapping technique. Trays and cases with lids may also be placed in an approved sterilization container with a sterilization sealing lid. Follow the sterilization container manufacturer's instructions for inserting and replacing sterilization filters in sterilization containers.

Instrument transport and storage cases and trays should be used under the following conditions:

- If possible, all devices should be disassembled before being placed in the case.
- All devices should be arranged so as to ensure the steam comes into contact with all the instrument surfaces. Instruments should not be stacked or placed in close contact (use storage holders or silicone mats provided).
- The user should check that there has been no spillage of the instrument case contents once the devices have been arranged in the case.
- Only devices manufactured and/or distributed by FH should be placed in the FH instrument tray.

Sterilization

- Reusable surgical instruments must always be disinfected according to the specified procedure prior to sterilization.
- The hospital is responsible for in-house instrument reassembly, inspection and packaging procedures following thorough cleaning so as to ensure sterilising steam penetration and suitable drying. Protection measures to be taken against sharp or potentially dangerous instrument edges should also be recommended by the hospital.
- Steam/moist heat sterilization is the preferred method recommended for FH orthopaedic instrument sets.
- The steriliser manufacturer's recommendations must **always** be followed. When sterilising several instrument sets in a sterilization cycle, take care not to exceed the maximum load specified by the manufacturer.
- The instrument sets should be prepared and packaged correctly in trays and/or cases to enable the steam to penetrate and come into direct contact with all the surfaces.
- Refer to the table below for the cycle sterilization settings validated by FH to provide a sterility assurance level (SAL) of 10⁻⁶.
- Ethylene oxide or plasma sterilization methods **must not** be used.

Validated steam sterilization settings:

| | Cycle type | Temperature (°C) | Temperature (°F) | Exposure time | Drying time |
|---|--------------|------------------|------------------|---------------|-------------|
| A | Prior vacuum | 132°C | 269.6°F | 4 minutes | 30 minutes |
| B | Prior vacuum | 134°C | 273.2°F | 18 minutes | 30 minutes |
| C | Prior vacuum | 134°C | 273.2°F | 3 minutes | 30 minutes |

Cycles B and C are not recommended in the United States, nor is it for the inactivation of prions.

This procedure was validated by an independent laboratory in accordance with AAMI TIR12.

Storage

- When handling packaged cases, take care to avoid damaging the sterile barrier. The healthcare establishment should define a storage shelf-life for packaged instruments, on the basis of the sterile packaging and the recommendations of the manufacturer of this packaging.

Responsibilities of the hospital with regards to FH instruments on loan / Number of potential uses

- FH Industrie does not claim a number of use of the instruments. Orthopaedic surgery instruments, if handled and used properly, can be used for several years. However, this service life may rapidly be reduced due to poor handling or unsuitable protection. Instruments that are no longer up to standards due to a long period of use, poor handling or unsuitable care should be returned to FH. Notify your FH representative of any problems with your instruments
- Before being returned to FH, FH instruments are required to undergo all the decontamination, cleaning, disinfection, inspection and final sterilization steps. The decontamination documentation should be furnished with the instruments returned to FH.

IMPORTANT INFORMATION FOR HOSPITAL PERSONNEL

Any serious accident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user has residence.

References

- AAMI TIR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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Referanslar

• AAMI TIR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.

• ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

• ISO 17665, Sterilization of health care products – Moist heat.

• ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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AAMI TIR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.

- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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Referenser.

- AAMI TIR12, Designing, testing, and labeling reusable medical devices for re processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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FR, FH ORTHO SAS