PROTHESES TOTALES DE GENOU TOTAL KNEE PROSTHESES TOTALE KNIETPROTHESEN PRÓTESIS TOTALES DE RODILLA PROTESI TOTALI DI GINOCCHIO PRÓTESES TOTAIS DO JOELHO VOLLEDIGE KNIETPROTHESEN PROTEZE TOTALE DE GENUNCH CAŁKOWITE ENDOPROTEZY STAWU

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en - Instructions for use TOTAL KNEE PROSTHESES
<p>WARNINGS AND INDICATIONS FOR USE PLEASE READ THIS DOCUMENT CAREFULLY</p> <p>Please request help from the distributor's sales department if you need additional information such as the surgical technique instructions which show details about implantation. In any case, refer to the documents provided.</p>

DESTINATION

"The implants mentioned in these instructions are intended to replace the knee joint during a total knee replacement."

DESCRIPTION

Preamble:

- A left TKR implant must be implanted on the left knee.
- A right TKR implant must be implanted on the right knee.
- A TKR implant that is not specified as being left- or right-sided may be implanted indiscriminately on a left or right knee.

For FHK TKR:

FH ORTHO offers a range of FHK total knee prostheses:

- fixed bearing (PF); posterior-stabilised (PS) or cruciate-retaining (PC).
- rotating platform mobile bearing (PR or "mobile").

The FHK left and right femoral components, sizes 1 through 7, are available in posterior-stabilised (PS) or standard versions, cemented (AC) or cementless (SC) versions.

The FHK tibial baseplates, sizes 1 through 7, are available in fixed (PF) or rotating bearing (PR or "mobile"), cemented (AC) or cementless (SC) versions.

FHK® tibial baseplates, FASTRACK version sizes 1 through 7, are available in fixed (PF) or rotating bearing (PR or "mobile" versions), cementless (SC)

The FHK inserts, sizes 1 through 7, thickness of 10, 12, 14, and 16 mm, are available (18 and 20 mm on special order) in the fixed bearing posterior-stabilised (PS) or cruciate-retaining (PC) version, or in the rotating version (PR or "mobile"), ultra-congruent (UC) or standard.

FHK TIBIA FULL PE TIBIAL COMPONENTS are cemented (AC) implants that can replace a baseplate and the corresponding PF insert. They are available in sizes 1 through 7, thickness 10, 12 ad 14 mm, in fixed bearing, cruciate-retaining (PC) or posterior-stabilised (PS) versions.

For FHK ASYMETRIC TKR:

FH ORTHO offers a range of FHK ASYMETRIC total knee prostheses, in the rotating version (PR or "mobile").

The left and right FHK ASYMETRIC tibial baseplates, sizes 1 through 6, are available in cemented (AC) or cementless (SC) versions.

The FHK ASYMETRIC inserts, sizes 1 through 6, thickness of 10, 12, 14, and 16 mm, are available (and 18 mm on special order).

FHK components are assembled with FHK ASYMETRIC baseplates and inserts to create the FHK ASYMETRIC TKR: details on compatibility between implants in the specific table for FHK ASYMETRIC TKR.

Implants shared by FHK TKR and FHK ASYMETRIC TKR:

The FHK patellas are available in a cemented (AC) resurfacing version or inlay version in the following dimensions:

FHK CEMENTED (AC) resurfacing patellas	CEMENTLESS (SC) FHK inlay patellas
Dimensions	Dimensions
Ø 30 / 8	Ø 22 / 7
Ø 34 / 8	Ø 25 / 7
Ø 34 / 10	
Ø 38 / 10	

The FHK diaphyseal stem extensions are available in the following dimensions:

FHK diaphyseal stem extensions
Dimensions
Ø 10 / L 70
Ø 14 / L 70
Ø 10 / L 110
Ø 14 / L 110

COMPATIBILITY BETWEEN DEVICES

The manufacturer does not authorise its product ranges to be combined with components from other manufacturers.

All the femoral components of our knee product ranges are lateralized for the left or right side.

Product combining is not possible between:

- FHK and FHK ASYMETRIC Rotating platform mobile bearing version (PR or "mobile"): between the size of the femoral component and the size of the rotating insert
- Fixed bearing version (PF): between the size of the fixed baseplate and the size of the PC or PS fixed bearing insert (see below)

The femoral component for the rotating insert is the same as the femoral component for the cruciate-retaining (PC) fixed bearing (PF) insert or for the TIBIAL COMPONENT / CRUCIATE-RETAINING / CEMENTED FHK TIBIA FULL PE. The posterior-stabilised (PS) femoral component is only compatible with the posterior-stabilised (PS) fixed bearing (PF) insert or the TIBIAL COMPONENT / POSTERIOR-STABILISED / CEMENTED FHK TIBIA FULL PE.

FHK TKR product range compatibility:

Materials	fixed bearing (PF) version					
	cruciate-retaining (PC) fixed bearing (PF)			posterior-stabilised (PS) fixed bearing (PF)		
	PF baseplate	PC PF insert	femoral component	PF baseplate	PS PF insert	PS femoral component
Compatibility between Sizes	1	1	1	1	1	1
			2			2
	2	2	1	2	2	1
			2			2
			3			3
	3	3	2	3	3	2
			3			3

4	4	4	4	4	4
		3			3
		4			4
5	5	5	5	5	5
		4			4
		5			5
6	6	6	6	6	6
		5			5
		6			6
7	7	7	7	7	7
		6			6

Materials	rotating platform mobile bearing version (PR or "mobile")		
	PR baseplate	PR or PR UC insert	femoral component
Compatibility between Sizes	1	1	1
	2	2	2
	3	3	3
	4	4	4
	5	5	5
	6	6	6
	7	7	7

FHK TKR with FHK TIBIA FULL PE TIBIAL COMPONENT:

Materials	FHK® TIBIA FULL PE TIBIAL COMPONENT			
	fixed bearing (PF) cruciate-retaining (PC)		fixed bearing (PF) posterior-stabilised (PS)	
	FHK® TIBIA FULL PE CRUCIATE-RETAINING TIBIAL COMPONENT	PC femoral component	FHK® TIBIA FULL PE POSTERIOR-STABILISED TIBIAL COMPONENT	PS femoral component
Compatibility between sizes	1	1	1	1
		2		2
	2	1	2	1
		2		2
		3		3
	3	2	3	2
		3		3
		4		4
	4	3	4	3
		4		4
	5	4	5	4
		5		5
		6		6
	6	5	6	5
		6		6
7	6	7	6	
	7		7	

FHK ASYMETRIC TKR product range compatibility:

- A left-sided FHK ASYMETRIC tibial baseplate must be assembled with a left-sided FHK ASYMETRIC insert
- A right-sided FHK ASYMETRIC tibial baseplate must be assembled with a right-sided FHK ASYMETRIC insert

Materials	FHK ASYMETRIC baseplate	FHK ASYMETRIC insert	FHK femoral component
	1	1	1
Compatibility between Sizes	2	2	2
	3	3	3
	4	4	4
	5	5	5
	6	6	6
	7	7	7

Combinations within the same product range are possible and more generally cemented and cementless implants are compatible with one another. The patellas (made of UHMWPE) are compatible with all sizes of femoral components. The diaphyseal stem extensions (made of CoCr) are compatible with all sizes of tibial baseplates.

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are noted on each product label, in particular the cementless (SC) coating materials. Stainless steel and chrome-cobalt alloys typically contain nickel. The implant is radiopaque unless it is made of UHMWPE.

PERFORMANCE

Total knee prostheses make it possible to reduce pain, restore the functional anatomy of the knee joint, and allow patients to return to activities that are compatible with a joint replacement.

INDICATIONS FOR USE:

Type of implant	FHK®	FHK® FASTRACK	FHK® FULL PE	FHK® ASYMETRIC
Painful and debilitating gonarthrosis caused by: osteoarthritis, rheumatoid arthritis, post-traumatic arthritis, and valgus, varus, and flexum type deformities	✓	✓	✓	✓
Aseptic necrosis of the femoral condyle	✓	✗	✗	✓

CONTRAINDICATIONS

- Infection, or latent infection
- A mental or neuromuscular disorder that might create an unacceptable risk of instability, prosthetic fixation failure, or complications after surgery
- Insufficient bone stock for proper implant fixation
- It is the surgeon's responsibility to make sure that the patient does not have any known allergies to one of the compounds of the material listed on the product label.
- Metabolic diseases that might compromise bone growth
- Drug addiction
- Incomplete bone growth

PRECAUTIONS

GENERAL PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.
- The outcome of arthroplasty depends on the patient's medical history (particularly orthopaedic antecedents). It is indispensable that the patient be psychologically prepared.
- Patients must be informed of the limitations of the prosthesis, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.
- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.
- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.
- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone and should remove excess cement.
- Care should be taken to ensure that the presence of other devices will not harm the functional integrity of the device.
- Do not combine implants from various origins or different systems provided by FH ORTHO. The only possible combinations of components are indicated above (§ Compatibility between devices)

CEMENTLESS IMPLANTS

During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants must not be cemented. The type of coating is identified on the label.

IMPLANTATION TECHNIQUE

Remember to comply with 4 imperatives:

- Respect the normal mechanical axis of the lower limb (Hip-Knee-Ankle angle of 180°)
- The stability of the joint obtained by carefully verifying ligament tension during flexion and extension should be carried out using the test devices provided for this purpose by the manufacturer. The surgeon will choose the appropriate implant accordingly after completing the tests.
- Good freedom of functional range of motion
- Precise restoration of the height of the joint space

The specific ancillary equipment may be delivered non-sterile, single-use or reusable, or sterile single-use (for example, saw blades). This equipment is necessary for implanting implants designed by FH industrie and distributed by FH ORTHO or its distributors. The surgical techniques are available from FH ORTHO or its distributors. For certain implants, FH ORTHO or its distributors may provide x-ray templates to enable the surgeon to evaluate the size of the implant.

CONDITIONS OF USE

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. The implantation must be performed with suitable non-damaged instrumentation by applying the recommendations in the surgical technique and instructions for use, to treat a patient with the indications defined above, either as a primary or revision arthroplasty procedure.

WARNINGS

Patients receiving a joint replacement must be informed that the lifespan of the implant may depend on their weight and their activity level. In obese patients, excess weight might put stress on the prosthesis and compromise either fixation of the device or the device; the use of a tibial pin (FHK diaphyseal stem extension) is particularly recommended in such cases (unless using a FULL PE tibial component).

- Damaged or defective implants must not be used.

- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics
- Using an implant that is of the wrong size or lateralisation runs the risk of decreasing its stress resistance.
- Components are part of a modular system. FH Industrie refuses any responsibility in the event that components from another manufacturer are used, or that its components are not used as they are intended to be.
- The implants must not be modified or subjected to any treatment.
- Restrictions of use: comply with the combinations authorised by the manufacturer

ADVERSE EFFECTS

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

- pain
- loss of mobility
- wear of the polyethylene components. Wear particles, especially from components made of UHMWPE, may cause osteolysis that may require another surgery.
- implant rupture as a result of inappropriate activities, trauma, or other stress specific to the patient's activities.
- loosening or dislocation of an implant following inadequate initial fixation, latent infection, early or excessive stress, poor positioning of the components, or trauma.
- allergy to one of the compounds of the material noted on the product label.
- possible evolving, sometimes asymptomatic, bone resorption that may appear around prosthetic components as a result of reactions to a foreign body.
- Hematoma, slowed healing, deep vein thrombosis, pulmonary thrombosis, vascular lesions
- Pulmonary embolism
- Infection, and possibly death
- A cardiovascular problem

These adverse effects may lead to another operation or a revision

IMPORTANT INFORMATION TO GIVE THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

Metal implants can move and/or heat up during an MRI exam. Metal implants can cause artifacts (signal distortion and/or signal loss) which may interfere with the interpretation of the imaging results. During MRI exams, the patient must inform the medical staff that he/she has a metallic implant.

PACKAGING AND STERILIZATION

The information on the product label ensures manufacturing traceability.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly closed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned.

The product must not be resterilised by any method whatsoever at the risk of altering its mechanical characteristics. The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation or green for ETO gas sterilisation. This colour may be altered by poor storage conditions: heat, humidity, light.

In all cases, an indicator showing one of the pre-sterilisation colours, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilised implants cannot be returned. The expiry date is noted on the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

The inserts, patellas and FULL PE tibial components are sterilised by:

The other components are sterilised by:

Date of first CE marking for FHK PF and PR 14/05/2009
 Date of first CE marking for FHK PR UC platform and inlay patella Ø22 16/01/2012
 Date of first CE marking for FHK T7 : 04/05/2016
 Date of first CE marking for FHK ASYMETRIC: 26/04/2016
 Date of first CE marking for FHK® tibial baseplate FASTRACK version: 14/01/2020
 Date of first CE marking for FHK® TIBIA FULL PE: 16/10/2020

The reference text is the French text.

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