FH ORTHOPEDICS
quality for health

FH ORTHOPEDICS S.A.S
3 rue de la Forêt - F 68990 HEIMSBRUNN
Tel.: +33 3 89 81 90 92 / Fax: +33 3 89 81 80 11
e-mail: orthopedie@fthorthopedics.fr
www.fthorthopedics.com

POLSKA, IMPLANTS INDUSTRIE
Ul. Garbary 95/A6,
61-757 Poznan
Tel.: +48 61 863 81 27 / Fax: +48 61 863 81 28
Email: fh.orthopedics@poczta.internetdsl.pl

USA, FH ORTHOPEDICS INC.
4118 N. Nashville Ave.
Chicago - IL 60634
Tel.: +1 (773) 290 1039 / Fax: +1 (773) 539 9328
e-mail: info-us@fthorthopedics.com

FH ORTHOPEDICS

ANCORAGE LIGAMENTAIRE CorTTape®
CorTTape® LIGAMENT ANCHORING
CorTTape® BANDVERANKERUNG
ANCLAJE LIGAMENTARIO CorTTape®
SISTEMA DI ANCORAGGIO PER LEGAMENTI
CorTTape®
ANCORAGEM LIGAMENTAR CorTTape®
LIGAMENTAIRE VERANKERING CorTTape®
MOCOWANIE WIĘZADŁA CorTTape®

CE 0459

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en

Instructions for use
CorTTape® LIGAMENT ANCHORING

WARNINGS AND INDICATIONS FOR USE PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical techniques which show details about implantation. In all cases, please consult the included documents.

DESTINATION

For implantation in the femoral section during ligament surgery on the knee.

DESCRIPTION

FH ORTHOPEDICS proposes a range of ligament anchoring systems for the knee. They must be used according to the surgical technique only.

The CorTTape® ligament anchoring range is available in two versions:
CorTTape® TENDON FIXATIONS: including the following elements:

- Titanium Plates
- Adjustable Loop Tape

CorTTape® SHIFTED TENDON FIXATIONS: including the following elements:

- Titanium Shifted Plates
- Adjustable Loop Tape

CorrTape® ligament anchoring

CorrTape® ligament anchoring	
CorrTape® TENDON FIXATIONS	✓
CorrTape® SHIFTED TENDON FIXATIONS	✓

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are indicated on each product label. Stainless steel and cobalt-chrome (CoCr) alloys typically contain nickel. The implant is radiopaque unless made from PEEK OPTIMA®, PET or Poly-L-lactide

COMPATIBILITY BETWEEN DEVICES

The ligament anchoring systems TLS® Ti, TLS® PEEK, TLS® BIO-C are designed for use with the TLS® tendon anchor tape only and are not therefore compatible with other systems, **except when used with a CorrTape® ligament anchoring system in the femur.**

TLS® Ti screw				
Diameter/Length	20	25	30	35
Ø 9	✓	✓	✓	✓
Ø 10	✓	✓	✓	✓
Ø 11	✓	✓	✓	✓
Ø 12	✓	x	x	x

TLS® Peek screw		
Diameter/Length	20	25
Ø 10	✓	✓

TLS® Bio-C screw		
Diameter/Length	20	25
Ø 10	✓	✓

The ancillary device to be used to put the implants from the ligament anchoring range into place is specified in the surgical technique.

PERFORMANCE

The CorrTape® ligament anchoring range secures the graft to the femur and TLS® ligament anchoring systems must be used on the tibia for reconstruction of the cruciate ligament in the knee.

INDICATIONS FOR USE:

Partial or total reconstruction of the anterior and/or posterior cruciate ligaments.

CONTRAINDICATIONS

- 1 - Infection, or latent infection.
- 2 - A mental or neuromuscular disorder that might create an unacceptable risk of instability, prosthetic fixation failure, or complications after surgery.
- 3 - Insufficient bone stock.
- 4 - Known allergy to any compound of the material mentioned on the product label.
- 5 - Metabolic diseases that might compromise bone regrowth.
- 6 - Uncooperative patient unable to follow recommendations.

PRECAUTIONS

The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use. The outcome of surgery depends on the patient's orthopaedic antecedents. It is indispensable that the patient be psychologically prepared.

Schroef TLS® Peek		
Diameter/ Lengte	20	25
Ø 10	✓	✓

Schroef TLS® Bio-C		
Diameter/ Lengte	20	25
Ø 10	✓	✓

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Patients must be informed of the limitations of the prosthesis, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised about how to modify their activities accordingly.

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.

- The components and all polished or coated surfaces must be protected against abrasion, scratches and any other damage by metal or abrasive objects.

- Care must be taken to ensure that the presence of other devices is no obstacle to the functional integrity of this system.

- Do not combine implants from various origins. The only possible combinations of components are indicated in the section on

IMPLANTATION TECHNIQUE

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FH ORTHOPEDICS or its distributors.

CONDITIONS OF USE

The ligament anchoring systems must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents the indications defined above.

WARNINGS

- Patients receiving an implant must be informed that the lifespan of the implant may depend on their weight and their level of activity.

- Damaged or defective implants must not be used.

- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.

- For the components of a modular system, FH Industrie refuses any responsibility in the event that components from another manufacturer are used, or that its components are not used as they are intended to be.

- The implants must not be modified or subjected to any treatment.

ADVERSE EFFECTS

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

- the risk of implant rupture as a result of inappropriate activities, trauma, or other stress specific to the patient's activities.

- the risk of an implant loosening following inadequate initial fixation, latent infection, early or excessive stress, poor positioning of the components, or trauma,

- the risk of allergy to one of the compounds of the material noted on the product label,

- side effects that may lead to another operation or a revision,

- possible evolving, sometimes asymptomatic, bone resorption that may appear around implants as a result of reactions to a foreign body,

- pain.

IMPORTANT INFORMATION TO GIVE THE SURGEON AND THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

Metallic implants carry a risk of migration and/or heating during an MRI exam. During MRI exams, the patient must inform the medical staff that he/she has a metallic implant.

PACKAGING AND STERILIZATION

The information on the product label ensures manufacturing traceability.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly sealed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be resterilised by any method whatsoever at the risk of altering its mechanical performance. The sterilization indicator sticker on the outer packaging should be red for radiation sterilization (R) or green for EIO gas sterilization (EO). This colour may be altered by poor storage conditions: heat, humidity, light etc. In all cases, an indicator showing one of the pre-sterilization colours, yellow or purple, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilised implants cannot be returned. The expiry date is noted on the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly ensuring environmental preservation and in strict compliance with sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING

CorrTape® TENDON FIXATIONS: 2015

CorrTape® SHIFTED TENDON FIXATIONS : 2015

The reference text is the French text.

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