



Technical sheet
REVISION OR ABLATION


ARROW I

GROUPE
FH ORTHO[™]

CASE n°1


Replacement of ARROW I
intermediate implants (identical)

Keep the stem
in place
ARROW I



The **ARROW I revision** ancillary **A266119** allows replacement
of intermediate parts (insert and glenosphere)

Keep the baseplate or
cemented glenoid **in place**
ARROW I



TO ORDER



ARROW II ancillary + **ARROW I Revision Ablation** ancillary
and associated implant nomenclatures



* order universal anatomical
ancillary ARROW I + the associated
implant nomenclature

CASE n°2

Change the stem and humeral intermediate implants
(humeral insert or humeral head) **in ARROW II**



ARROW II



Keep the baseplate or
cemented glenoid **in place**

ARROW I



Replacement of glenoid
intermediate implants
(insert or glenosphere) in

ARROW I



TO ORDER




ARROW II ancillary + **ARROW I Revision Ablation** ancillary
and associated implant nomenclatures

* order universal anatomical
ancillary ARROW I + the associated
implant nomenclature


CASE n°3

Change of baseplate and glenoid intermediate implants (glenosphere or glenoid insert) **in ARROW II**



Keep the stem
in place
ARROW I



Replacement of humeral
intermediate implant
(humeral insert or humeral head) in
ARROW I



ARROW II




TO ORDER



ARROW II ancillary + **ARROW I Revision Ablation** ancillary
and associated implant nomenclatures

CASE n°4

Change the whole platform in **ARROW II**




TO ORDER



ARROW II ancillary + **ARROW I Revision Ablation** ancillary
and associated implant nomenclatures

ARROW I IMPLANTS



ARROW I REVISION ABLATION ANCILLARY (ref. A 266 119)

ARROW I Revision implant nomenclature (ref. I 272 198)



ARROW I HUMERAL INSERT

Ref.	Designation
257 057	ARROW humeral insert D36 H00
257 058	ARROW humeral insert D36 H05
257 059	ARROW humeral insert D36 H10
257 060	ARROW humeral insert D39 H00
257 061	ARROW humeral insert D39 H05
257 062	ARROW humeral insert D39 H10
257 063	ARROW humeral insert D42 H00
257 064	ARROW humeral insert D42 H05
257 065	ARROW humeral insert D42 H10



ARROW I HUMERAL HEAD

Ref.	Designation
260 532	ARROW centred humeral head D48 H18
260 535	ARROW centred humeral head D50 H19
260 550	ARROW centred humeral head D54 H19




ARROW I GLENOSPHERE

Ref.	Designation
260 566	ARROW glenosphere D36
260 567	ARROW glenosphere D39
260 568	ARROW glenosphere D42

ANATOMICAL METAL-BACK PROSTHESIS REVISION ARROW I

ARROW I ANATOMICAL UNIVERSAL BASKET (ref. A 272 069)

ARROW I implant nomenclature Anatomical revision (ref. I 272 104)



CENTRED HUMERAL HEAD

Ref.	Diameter	Height (mm)
260 537	40	15
260 538	40	17
260 539	44	16
260 540	44	18
260 541	46	16
260 542	46	18
260 543	46	21
260 544	48	16
260 545	48	18
260 546	48	21
260 547	50	17
260 548	50	19
260 549	50	21
260 550	54	19
260 551	54	21



OFF-CENTRED HUMERAL HEAD

Ref.	Diameter	Height (mm)
260 526	44	16
260 527	44	18
260 528	46	16
260 529	46	18
260 530	46	21
260 531	48	16
260 532	48	18
260 533	48	21
260 534	50	17
260 535	50	19
260 536	50	21



GLENOID INSERT

Ref.	Size
260 556	44
260 557	46
260 558	48

GROUPE
FH ORTHOTM



DISTRIBUTEURS
DISTRIBUTORS

FR, FH ORTHO SAS
3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHOPEDICS INC.
OrthoEx
7327 E Tierra Buena Lane
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fhortho-us.com
www.fhortho.com

PL, FH ORTHO POLSKA
Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com



FABRICANT
MANUFACTURER

FR, FH INDUSTRIE
6 rue Nobel, Z.I. de Kernevez
29000 QUIMPER - FRANCE
Tél. +33 (0)2 98 55 68 95
Fax : +33 (0)2 98 53 42 13
contact-fhi@fhortho.com
www.fhortho.com