



45b screw



Bevelled head

Minimisation of screw protrusion and soft tissue irritation thanks to the bevelled profile of the screw head



Threading

Entirely threaded, non-compressive



Power driven screw

Self-tapping
Self-drilling



Specific instrumentation

Screwdriver with bevel identifier



Complete program

Available from length 16 to 60 mm

INDICATIONS

- Treatment of hallux valgus of the 1st ray with osteosynthesis without compression,
- Osteosynthesis without compression of forefoot fractures.

IMPLANT REFERENCES

Ref.	Designation
270 681	bePOD bevelled cannulated 45B screw Ø3 LG16 mm
270 682	bePOD bevelled cannulated 45B screw Ø3 LG18 mm
270 683	bePOD bevelled cannulated 45B screw Ø3 LG20 mm
270 684	bePOD bevelled cannulated 45B screw Ø3 LG22 mm
270 685	bePOD bevelled cannulated 45B screw Ø3 LG24 mm
270 686	bePOD bevelled cannulated 45B screw Ø3 LG26 mm
270 687	bePOD bevelled cannulated 45B screw Ø3 LG28 mm
270 688	bePOD bevelled cannulated 45B screw Ø3 LG30 mm
267 823	bePOD bevelled cannulated 45B screw Ø3 LG32 mm
267 824	bePOD bevelled cannulated 45B screw Ø3 LG34 mm
267 825	bePOD bevelled cannulated 45B screw Ø3 LG36 mm
267 826	bePOD bevelled cannulated 45B screw Ø3 LG38 mm

Ref.	Designation
267 827	bePOD bevelled cannulated 45B screw Ø3 LG40 mm
267 828	bePOD bevelled cannulated 45B screw Ø3 LG42 mm
267 829	bePOD bevelled cannulated 45B screw Ø3 LG44 mm
267 830	bePOD bevelled cannulated 45B screw Ø3 LG46 mm
267 831	bePOD bevelled cannulated 45B screw Ø3 LG48 mm
267 832	bePOD bevelled cannulated 45B screw Ø3 LG50 mm
270 689	bePOD bevelled cannulated 45B screw Ø3 LG52 mm
270 690	bePOD bevelled cannulated 45B screw Ø3 LG54 mm
270 691	bePOD bevelled cannulated 45B screw Ø3 LG56 mm
270 692	bePOD bevelled cannulated 45B screw Ø3 LG58 mm
270 693	bePOD bevelled cannulated 45B screw Ø3 LG60 mm

INSTRUMENTATION




267 807 - Motorized screwdriver for 45B screw



267 808 - Manual screwdriver for 45B screw



267 809 - Drill for 45B screw



267 811 - Measurement tool for 45B screw

254 920 - K-wire for cannulated screw Ø 1 mm lg 130 mm

269 057 - Transmetatarsal K-wire Ø 1.5 mm lg 150 mm



Manufacturer: FH INDUSTRIE

6 rue Nobel 29000 Quimper - France
e-mail : contact-fhi@fhortho.com

GROUP
FH
ORTHO

FR, FH ORTHO SAS

3 rue de la Forêt - Zone Industrielle
BP 50009
68990 HEIMSBRUNN CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHO INC

4908 N. Elston
Chicago - IL 60630 - USA
Phone: +1 (773) 290 1039 / 844-77 FHINC
Fax: +1 (708) 667 7618
info-us@fhorthopedics.com
www.fhortho.com

UK, FH ORTHO LTD

Unit 1b, Century Park, Valley way
Swansea Enterprise Park
Swansea, SA6 8RP - UK
Phone: +44 (0) 1792 464792
Fax: +44 (0) 844 412 7674
customer-servicesUK@fhorthopedics.com
www.fhortho.com

PL, FH ORTHO POLSKA

Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com