

FOOT

tictac

TOE INTERPHALANGEAL CAGE
TOE ARTHRODESIS CAGE




TICTAC Ø2,5 - 0°
ref. 270 390



TICTAC Ø2,8 - 0°
ref. 268 612



TICTAC Ø3,1 - 0°
ref. 268 614



TICTAC Ø3,4 - 0°
ref. 268 616



TICTAC Ø2,5 - 10°
ref. 270 391



TICTAC Ø2,8 - 10°
ref. 268 613



TICTAC Ø3,1 - 10°
ref. 268 615



TICTAC Ø3,4 - 10°
ref. 268 617



TICTAC Ø2,5 - 20°
ref. 270 392



TICTAC Ø2,8 - 20°
ref. 270 393



TICTAC Ø3,1 - 20°
ref. 270 394



TICTAC Ø3,4 - 20°
ref. 270 395

- > **Designed to fit all the patient's anatomy:**
 - 3 angulations: 0°/10°/20°
 - 4 diameters: 2,5mm, 2,8mm, 3,1mm and 3,4mm

- > **Intra-medullar press-fit** achieved by the implant design.

- > **Made of PEEK** as a way to improve the biocompatibility of load-bearing implants.

- > **Dedicated instruments**

> Instrumentation

DRILL Ø2

ref. 266 818



IPP CAGE TICTAC SHAPER P1 Ø2,5 BEPOD

ref. 271 066



IPP CAGE TICTAC SHAPER P2 Ø2,5 BEPOD

ref. 271 067



IPP CAGE TICTAC SHAPER P1 Ø2,8 BEPOD

ref. 268 618



IPP CAGE TICTAC SHAPER P2 Ø2,8 BEPOD

ref. 268 619



IPP CAGE TICTAC SHAPER P1 Ø3,1 BEPOD

ref. 268 620



IPP CAGE TICTAC SHAPER P2 Ø3,1 BEPOD

ref. 268 621



IPP CAGE TICTAC SHAPER P1 Ø3,4 BEPOD

ref. 268 622



IPP CAGE TICTAC SHAPER P2 Ø3,4 BEPOD

ref. 268 623




IPP CAGE TICTAC DRIVER HOLDER BEPOD

ref. 268 624



HOLDER BEPOD

ref. 271 147



FR, FH ORTHO SAS

3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fortho.com
www.fortho.com

USA, FH ORTHO INC.

OrthoEx
16000 North 80th Street
Suite E
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fortho-us.com
www.fortho.com

UK, FH ORTHO LTD

Unit 1b, Century Park, Valley way
Swansea Enterprise Park
Swansea, SA6 8RP - UK
Phone: +44 (0) 1792 464792
Fax: +44 (0) 844 412 7674
customer-servicesUK@fortho.com
www.fortho.com

PL, FH ORTHO POLSKA

Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fortho.com

FR, FH INDUSTRIE

6 rue Nobel, Z.I. de Kernevez
29000 QUIMPER - FRANCE
Tél. +33 (0)2 98 55 68 95
Fax : +33 (0)2 98 53 42 13
contact-fhi@fortho.com
www.fortho.com

